# MED-RX® Oral/Enteral Syringe with ENFit® Connector, Non-sterile, Instructions for Use

## INDICATIONS

This device is designed for use as a dispenser, a measuring device, and a fluid transfer device. It is used to deliver fluids to patients through enteral extension sets and feeding tubes or orally. This product is designed to connect with ENFit® connectors. It is intended for use by medical professionals on all patient populations.

## PRODUCT CODES

Product Code	Description
54-4010	10CC ORAL SYRINGE NON-STERILE, ENFit® ENTERAL SYRINGE
54-4020	20CC ORAL SYRINGE NON-STERILE, ENFit® ENTERAL SYRINGE
54-4030	30CC ORAL SYRINGE NON-STERILE, ENFit® ENTERAL SYRINGE
54-4060	60CC ORAL SYRINGE NON-STERILE, ENFit® ENTERAL SYRINGE

Table 1

## CONTRAINDICATIONS

- The syringe is not intended for parenteral use.
- The syringe is not designed for use with Luer syringes, including Luer Lock and Luer Slip, or friction fit connections.

## CAUTIONS

- Oral/enteral syringes are for one time, use only.
- Tip caps are a choking hazard. Keep out of reach of children.

## DIRECTIONS:

- Use Aseptic Technique.
- Follow medical professional instructions or hospital protocol in opening package and use of contents.
- Fully depress the plunger before filling oral/enteral syringe. Draw liquid (for feeding) into syringe.
- To expel excess air bubbles, hold syringe with tip pointing up. Flick the side of syringe to cause bubbles to rise toward tip. Express the bubbles into the medication vial or straw. Repeat until bubbles are eliminated.
- Disconnect the syringe from the filling connector and place the orange cap (packaged separately) on the syringe and label until ready to use (label not included).
- Administration should only be accomplished using ENFit® sets and feeding tubes with a ENFit® male connector.

Single Use Only. Do not reuse.

Lot Information

Non-sterile

Not made with natural rubber latex.

Not made with di (2-ethyl hexyl) phthalate (DEHP).

Not made with bisphenol A (BPA).

ENFit® is a federally registered trademark of GEDSA, used with permission.

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