



TAMPER EVIDENT CAP FOR ENFIT® SYRINGES

The Prep-Lock™ Tamper Evident Cap for ENFit® Syringes is specifically designed to meet the ISO 80369-3 Standard for ENFit Connectors. It is also the only ENFit Syringe Cap available with built-in tamper-evident technology.

Compatible with a low dose and standard ENFit Syringes, they provide an additional level of security to syringe contents during transportation, prior to administration.

Once the Tamper Evident Cap is installed, it cannot be removed without evidence of access and cannot be reinstalled after removal.

- **High Compatibility**

Compatible with both standard and low dose ENFit™ syringes.

- **Actively Deter Diversion**

Tampering is immediately evident when outer sleeve is detached or missing

- **Enhance Drug Security**

Elevate your standard of care with one product to guard all your ENFit® Syringes.

- **Ensure Integrity**

Help ensure the accuracy and integrity of syringe contents.

Features:

Packaged non sterile.

This product does not contain DEHP nor BPA.

This product is not made with natural rubber latex.

Made with Polypropylene & Polystyrene.

ORDER #	COLOR	QTY
57-220	PURPLE	BOX OF 1000
57-220-100	PURPLE	BOX OF 100

100 UNITS PER BAG, SOLID OUTER SLEEVE

ENFit® is a GEDSA trademark.